

Sound Wave - Personal

Anywhere at any time~

You can enjoy **Acoustic vibration** massager or cleansing with
Sound wave – Personal

Features

- Anywhere at anytime you can enjoy Acoustic vibration massage
- Straight for ward with 3 modes and intensity control button
- Available to massage from crown to toe in one device with 7 types of massage handle
- Facial cleansing and scalp cleansing are available – Waterproof

How to use

- Apply cleansing pad or appropriate massage handle and select the features that you want.
- Then gently massage or cleanse your face, scalp or body



Sound Wave - Personal

Anywhere at any time~
You can enjoy **Acoustic vibration** massager or cleansing with
Sound wave – Personal



Relieve dolorific point

Accelerate Nerve cells & Metabolism

Musculo-skeletal disorders mitigation

Relief cold hands and feet & face flush

Sliming face & body

Scalp cleansing with shampoo

Increase the circulation of the lymph system & Edema relief



Sound Wave - Personal

Anywhere at any time~
You can enjoy **Acoustic vibration** massager or cleansing with
Sound wave – Personal



Product contents

<u>Name</u>	Sound wave personal	<u>Power supply unit</u>	AC 100~240V
<u>Model</u>	NDS-H100		DC 5V/1A Cable length 1.5M
<u>Weight</u>	1.1kg (Body part: 195g)	<u>Charger</u>	Lithium-ion battery 500mAh Charging time 90Min / Using time 60 M
<u>Size</u>	<u>body part</u> 73x130x70	<u>Display</u>	LED
	<u>mm rests</u> 85x110x45		
	mm		
Authentication number (KC) : R-REM-NFA-NDS-H100			




Sound Wave - Personal

Anywhere at any time~
You can enjoy **Acoustic vibration** massager or cleansing with
Sound wave – Personal




SCALP MASSAGE – Scalp wire head



- ① Repeat the procedure of pressing the front of the head
- ② Repeat the process of pressing the top of the head
- ③ Repeat the procedure for pressing the back of the head
- ④ Repeat the procedure for pressing the right and left side of the head in regular order

Helpful for the scalp health by promoting the blood circulation and prevention of hair loss (10 minutes for each parts)


SCALP(HAIR) SHAMPOO OR TREATMENT MASSAGE – Scalp silicon head



- ① Press your hair as if it were a comb from the forehead to the top of your head and massage and massage
- ② Massage from the neck to the top of the head with a comb from the bottom to the back of the neck
- ③ Apply the pressure from the border of the forehead to the crown and back of the head
- ④ At the center of the line, press and comb the hair and massage the front, side and back of the head.
- ⑤ Slowly massage like a comb according to hair

Helpful for the scalp and hair health by deep cleansing while shampooing.

FACIAL MASSAGE – Disc head



Message time : 5~ 10 minute
(MESSAGE WITH AMPLIFIER ESSENCE)

- ① **FROM INSIDE TO OUTSIDE**
Drag the facial line from the center of the chin to the bottom of the ear, Drag from the side of the nose to the middle of the ear, Drag from under your eyes to your temples
- ② **FROM BOTTOM TO TOP**
From the center of the chin, go through the lip line to the bottom of the eye, For the chin part, overall, lift the face and pull it up vertically
- ③ **SMOOTH OUT WRINKLES AROUND THE EYES**
Drag the area with wrinkles from the bottom to the top intensively
- ④ **LIFTING THE FOREHEAD**
Gradually move from the eyebrow up to hair line.


Helpful for the anti-aging just like smoothing out the wrinkles and provide elasticity to the skin.

Sound Wave - Personal

Anywhere at any time~
You can enjoy **Acoustic vibration** massager or cleansing with
Sound wave – Personal




FACIAL CLEANSING DISC – Head+Disc silicon



- 1 First attach the disc silicon to the disc head then connect to the head part of the body device.
- 2 Wash your face with , put the foam cleanser on the palm of your hand to make enough bubbles.
- 3 Apply a light layer of foam evenly around the eyes
- 4 Like massaging the foam applied, cleansing the end of chin and line of chin round 20 seconds with the disc silicon connected to the body device.
- 5 Cleansing the T-zone area and the area where has lot of sebum (forehead and nose) about 20 seconds.
- 6 Finally, cleansing each cheek for 10 seconds

Helpful for the pore and exfoliation by the effect of deep cleansing

NECK MASSAGE T-head (Myofascial) & Circular cone-Head (Acupoint, Blood point)



1 LYMPH CIRCULATION

Start with the ear's lymph node and lower it down along the neckline to the shoulder then go down to the neck again

2 VERTICAL CARE

Moves up and down from neck to chin vertically.

3 BACK NECK LINE ACUPRESSURE

Press deeply along the hair area and slowly apply acupressure. You can feel the lightness of your head through massaging the whole acupressure


4 MAKING A NECKLINE

Loosen the muscles by sweeping it from the back of the neck to the shoulder. Neckline becomes longer and makes a beautiful neckline.

Helpful for the elasticity of skin and relieve a ring on the skin called the growth of a tree

BODY MASSAGE POINTS

T-head (myofascial, lymph) →
circular cone head (acupoint, blood point) ●
ball-head (joint, ligament, tendon) ★



It is recommended to massage the desired area slowly according to lymph nodes in the direction of the arrow.
Recommended to massage the entire body from top to bottom over time. (Do not use over 10 minutes)